












INNOCOR® Disposables/Consumables for INNOCOR LCI

INNOCOR® Model INN00400-1b Disposables/Consumables		Picture
Part no.	Part description	
BAG00001	<p>Rebreathing Bag, 3 L – Latex free To be used with adaptor MEC00120 (included in the STARTKIT) Recommended use: 3 per year</p>	
BAG00008	<p>Rebreathing Bag, 2 L – Latex free To be used with adaptor MEC00120 (included in the STARTKIT) Recommended use: 3 per year</p>	
CLI00002	<p>Nose Clips Sold in multiples of 5 pieces only, the price is per item Recommended use: new per each patient</p>	
FIL00001	<p>Particle Filters for Gas Analyzer Sold in multiples of 5 pieces only, the price is per item Recommended use: change every 500 tests - at least once a year (part of recommended yearly calibration)</p>	
FIL00103	<p>FILTER: HYGROBOY ISO 22M/15F-22F/15M for subjects < 10 kg Recommended use: new per each patient</p>	
FIL00105	<p>FILTER: BARRIERBAC S - Electrostatic Filter for subjects > 10 kg Recommended use: new per each patient</p>	
FIL00201	<p>Nylon Mesh for Flowmeter - Pediatric Sold in multiples of 5 pieces only Recommended use: to be changed once a week/month depending on the number of measurements performed</p>	

<p>MOU00004</p>	<p>Mouthpieces 22 mm (Hans Rudolph, reusable) Sold in multiples of 2 pieces only, the price is per item Recommended use: reusable, simple cleaning instructions are provided, mild household cleaning solution</p>	
<p>TUB00102</p>	<p>Nafion Dryer Gas Tube 36" Recommended use: change every 500 test - at least once a year (part of recommended yearly calibration)</p>	
<p>GAS00001</p>	<p>Rebreathing Test Gas with Cylinder Gas cylinder contains approx. 150-200 tests depending on patient age and size</p>	
<p>SUB00064</p>	<p>Compact RVU Insert Recommended use: Silicone valve insert: change every 4 months or after a gas cylinder is used</p>	
<p>ABS00002</p>	<p>Disposable Scrubber (CO2 Absorber) Recommended use: Useable for approx. 20 manoeuvres</p>	